


PUBLICATION SPOTLIGHT

Sentimag® – Magseed®


Targeted Axillary Dissection (TAD) – results from 900+ target node localisations


Breast cancer patients diagnosed with clinically positive lymph nodes often receive neoadjuvant systemic therapy (NAST) with the aim to downstage or eradicate the disease. To accurately assess the response in the breast and the axilla, it is important that both the positive lymph node/s and the breast lesion are marked before NAST to be able to locate them later on. Systemic therapy can negatively impact lymphatic drainage and hence reduce the accuracy of the sentinel lymph node biopsy (SLNB). However, when SLNB is paired with removal of the previously positive target lymph node, a technique called Targeted Axillary Dissection (TAD), the operation becomes a lot more accurate. Studies of TAD have demonstrated a false negative rate of < 2% compared to SLNB alone. (Caudle *et al.* JCO 2016)

The Magseed® marker has already been used to safely and effectively localise over 100,000 breast lesions and axillary lymph nodes. Its small size makes it the ideal marker for TAD and can now be implanted in any soft tissue for as long as required. Clinical studies involving over 900 patients globally have demonstrated accurate placement, no migration once implanted and accurate removal. In addition to improving clinical outcomes, it also offers a better patient experience and promotes seamless radiology and OR scheduling. When the Magseed® marker is used in combination with the lymphatic tracer Magtrace®, they offer the world's only wire-free, radiation-free solution for TAD in one platform.

Placement success

**Longest placement
to surgery time**


Range = 0 – 172 days

Retrieval rate

Nodal identity

**Sentinel
node**

**Non-sentinel
node**

Analysis based on 19 publications

Clinical study results – selection

First author	Patients	Seeds placed	Placement success	Longest placement to surgery time	Retrieval rate	TLN = SLN
Martinez ¹	81	82	100 %	ND	100 %	81.5 %
Reitsamer ²	40	50	100 %	ND	100 %	65 %
Thill ³	12	14	100 %	172 days (average)	100 %	100 %
Salazar Gomez ⁴	28	29	100 %	ND	100 %	52 %
Greenwood ⁵	35	38	97 %	31 days	97 %	ND
Sinnett ⁶	15	15	93 %	ND	93 %	86 %
Miller ⁷	ND	129	100 %	ND	100 %	ND
Simons ⁸	50	50	100 %	30 days	100 %	80 %

Target lymph node localisation – clinical results

Publications

[1] Martínez et al. (2022): Evaluation of Axillary Lymph Node Marking with Magseed® before and after Neoadjuvant Systemic Therapy in Breast Cancer Patients: MAGNET Study. *The Breast Journal*. Volume 2022, ID 6111907. [open access]



[2] Reitsamer R et al. (2021): The applicability of Magseed® for Targeted Axillary Dissection in breast cancer patients treated with neoadjuvant chemotherapy. *The Breast* 57: 113–7. [abstract]

[3] Thill et al. (2020): Magseed®-guided long-term marking of target lymph nodes in neoadjuvant therapy of early breast cancer patients - first experiences and prospectives. *Senologie – Zeitschrift für Mammadiag. und -therapie* 17(02): e44–e45. [abstract]



[4] Salazar Gomez et al. (2019): Utilidad de la semilla magnética para la localización de los ganglios axilares tras tratamiento neoadjuvante. *4th Spanish Breast Congress*. [open access]

[5] Greenwood H et al. (2019): Feasibility of Magnetic Seeds for Pre-Operative Localization of Axillary Lymph Nodes. *AJR* 213: 1–5. [abstract]



[6] Sinnott V (2019): Magnetic seeds: An attractive localisation option for management of axillary node positive breast cancer. *EJSO* 45(5): 889. [abstract]

[7] Miller et al. (2019): Hospital system rollout and initial experience with stainless magnetized seeds for breast and lymph node localization. *Ann. Surg. Oncol.* 26(suppl. 1). [abstract]



[8] Simons JM et al. (2019): Prospective Trial of Magnetic Seed Localization of Clipped Nodes after Neoadjuvant Chemotherapy in Node Positive Breast Cancer. *Ann. Surg. Oncol.* 26 (suppl. 1). [open access]

Last update October 2022

Endomag®, Sentimag® and Magseed® are registered European Union trade marks of Endomag Ltd · www.endomag.com

Distributor EMEA: Sysmex Europe SE
Bornbarch 1, 22848 Norderstedt, Germany · Phone +49 40 52726-0 · Fax +49 40 52726-100 · oncology@sysmex-europe.com · www.sysmex-europe.com

Manufacturer: Endomag Ltd
330 Cambridge Science Park, Milton Road, Cambridge, CB4 0WN, United Kingdom

You will find your local Sysmex representative's address under www.sysmex-europe.com/contacts