


PUBLICATION SPOTLIGHT

Insights about a patient's immune response from a complete blood count

A blood count is part of almost every medical check-up and, as such, one of the most frequently ordered laboratory tests. With regards to infection diagnostics, the WBC count and differential accompanied by morphology information from a microscopy are a traditional source of early information about a patient's condition.

Advances in blood count technology bring new parameters and allow deeper insights into the blood cells than the cell counts only. Sysmex's XN-Series analysers offer a set of parameters (Extended Inflammation Parameters) that makes it possible to quantitatively assess immature granulocytes (IG) and the activation status of neutrophils (NEUT-RI, NEUT-GI) and lymphocytes (RE-LYMP, AS-LYMP). They are readily available from a routine complete blood count with differential.

Several research studies show how Extended Inflammation Parameters (EIP) can support infection diagnostics by measuring WBC immune cell activation from a CBC + DIFF.



List of references

Publications

[1] **Cetin N et al. (2023):** Immature granulocytes as biomarkers of inflammation in children with predialysis chronic kidney disease. *Pediatr Nephrol.* 38(1): 219–225.



[2] **Lemkus L et al. (2022):** The utility of extended differential parameters as a biomarker of bacteremia at a tertiary academic hospital in persons with and without HIV infection in South Africa. *PLoS One* 2022 Feb 17; 17(2).



[3] **Rutkowska E et al. (2021):** Usefulness of the New Hematological Parameter: Reactive Lymphocytes RE-LYMP with Flow Cytometry Markers of Inflammation in COVID-19. *Cells* 2021 Jan 6; 10(1): 82.



[4] **Rolla R et al. (2020):** Reduced activity of B lymphocytes, recognised by Sysmex XN-2000™ haematology analyser, predicts mortality in patients with coronavirus disease 2019. *Int J Lab Hematol* 2021 Feb; 43(1): e5–e8.



[5] **Linszen J et al. (2020):** A novel haemocytometric COVID-19 prognostic score developed and validated in an observational multicentre European hospital-based study. *Elife* 2020 Nov 26; 9: e63195.



[6] **Huang Y et al. (2019):** Immature granulocytes: A novel biomarker of acute respiratory distress syndrome in patients with acute pancreatitis. *J Crit Care* 50: 303–308.



[7] **Stiel L et al. (2019):** First visualization of circulating neutrophil extracellular traps using cell fluorescence during human septic shock-induced disseminated intravascular coagulation. *Thromb Res* 2019 Nov; 183: 153–158.



[8] **Henriot I et al. (2017):** New parameters on the hematology analyzer XN-10 (Sysmex TM) allow to distinguish childhood bacterial and viral infections. *Int J Lab Hematol* 2017 Feb; 39(1): 14–20.



[9] **Weimann K et al. (2015):** Intensive Care Infection Score-A new approach to distinguish between infectious and noninfectious processes in intensive care and medicosurgical patients. *J Int Med Res* 2015 Jun; 43(3): 435–51.



[10] **Nierhaus A et al. (2013):** Revisiting the white blood cell count: immature granulocytes count as a diagnostic marker to discriminate between SIRS and sepsis – a prospective, observational study. *BMC Immunol* 2013 Feb 12; 14:8.

