

**Trade name:** CELLCLEAN

**Current version :** 3.1.1, issued: 04.05.2022

**Replaced version:** 3.1.0, issued: 01.06.2021

**Region:** GB

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

### 1.1 Product identifier

**Trade name**

**CELLCLEAN**

### 1.2 Relevant identified uses of the substance or mixture and uses advised against

**Relevant identified uses of the substance or mixture**

For in-vitro diagnosis

**Uses advised against**

No data available.

### 1.3 Details of the supplier of the safety data sheet

**Address**

SYSMEX LTD  
Sysmex House  
Garamonde Drive  
Wymbush  
Milton Keynes  
Buckinghamshire  
United Kingdom

MK8 8DF

Telephone no. +44 (0)333 320 3460

Fax no. +44 (0)1908 267 901

**Advice on Safety Data Sheet**

info@sysmex-europe.com

**Identification of the manufacturer / supplier**

**Address**

SYSMEX EUROPE GMBH  
Bornbarch 1  
D-22848 Norderstedt

### 1.4 Emergency telephone number

For medical advice (in German and English):

+49 (0)551 192 40 (Giftinformationszentrum Nord)

## SECTION 2: Hazards identification

### 2.1 Classification of the substance or mixture

**Classification in accordance with Regulation (EC) No 1272/2008 (CLP)**

Aquatic Acute 1; H400

Aquatic Chronic 2; H411

Eye Dam. 1; H318

Skin Corr. 1; H314

**Classification information**

Product is classified as "Corrosive" based on the extreme pH-value, see:

- Regulation 1272/2008 (CLP), Annex. I, number 3.2.2.2 / 3.2.3.1.2

This product is assessed and classified using the methods and criteria below referred to in Article 9 of Regulation (EC) n° 1272/2008:

Physical hazards: determined through assessment data based on the methods or standards referred to in part 2 of Annex I to CLP

Health hazards and environmental hazards: determined through toxicological and ecotoxicological assessment data based on the methods or standards referred to in Part 3, 4 and 5 of Annex I to CLP.

### 2.2 Label elements

**Labelling according to Regulation (EC) No 1272/2008 (CLP Regulation)**

**Hazard pictograms**

Trade name: CELLCLEAN

Current version : 3.1.1, issued: 04.05.2022

Replaced version: 3.1.0, issued: 01.06.2021

Region: GB



GHS05



GHS09

**Signal word**

Danger

**Hazard statement(s)**

H314 Causes severe skin burns and eye damage.  
 H410 Very toxic to aquatic life with long lasting effects.

**Hazard statements (EU)**

EUH031 Contact with acids liberates toxic gas.

**Precautionary statement(s)**

P260 Do not breathe mist/vapours/spray.  
 P273 Avoid release to the environment.  
 P280 Wear protective gloves/protective clothing/eye protection/face protection.  
 P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water [or shower].  
 P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
 P310 Immediately call a POISON CENTER/doctor.  
 P391 Collect spillage.

**2.3 Other hazards**

PBT assessment  
 No data available.

vPvB assessment  
 No data available.

**SECTION 3: Composition/information on ingredients****3.1 Substances**

Not applicable. The product is not a substance.

**3.2 Mixtures****Hazardous ingredients**

No	Substance name	Additional information	
	CAS / EC / Index / REACH no	Classification (EC) 1272/2008 (CLP)	Concentration %
1	<b>sodium hypochlorite, solution</b>		
	7681-52-9 231-668-3 017-011-00-1 -	Aquatic Acute 1; H400 Aquatic Chronic 1; H410 EUH031 Eye Dam. 1; H318 Skin Corr. 1B; H314	>= 5.00 - < 10.00 wt%

Full Text for all H-phrases and EUH-phrases: pls. see section 16

No	Note	Specific concentration limits	M-factor (acute)	M-factor (chronic)
1	B	EUH031: C >= 5%	M = 10	M = 1

Full text for the notes: pls. see section 16 "Notes relating to the identification, classification and labelling of substances ((EC) No 1272/2008, Annex VI)".

**3.3 Other information**

5 % Cl active

**SECTION 4: First aid measures****4.1 Description of first aid measures****General information**

**Trade name:** CELLCLEAN**Current version :** 3.1.1, issued: 04.05.2022**Replaced version:** 3.1.0, issued: 01.06.2021**Region:** GB

In case of persisting adverse effects, consult a physician. Remove contaminated clothing and shoes immediately, and launder thoroughly before reusing.

**After inhalation**

Remove affected person from the immediate area. Ensure supply of fresh air.

**After skin contact**

Wash off immediately with soap and water.

**After eye contact**

Separate eyelids, wash the eyes thoroughly with water (15 min.). Call a doctor immediately.

**After ingestion**

Call a doctor immediately. Do not induce vomiting. Never give anything by mouth to an unconscious person.

**4.2 Most important symptoms and effects, both acute and delayed**

No data available.

**4.3 Indication of any immediate medical attention and special treatment needed**

No data available.

**SECTION 5: Firefighting measures****5.1 Extinguishing media****Suitable extinguishing media**

Product itself is non-combustible; adapt fire extinguishing measures to surrounding areas.

**Unsuitable extinguishing media**

No data available.

**5.2 Special hazards arising from the substance or mixture**

In the event of fire, the following can be released: Carbon dioxide (CO<sub>2</sub>); Carbon monoxide (CO); Chlorine (Cl<sub>2</sub>)

**5.3 Advice for firefighters**

Use self-contained breathing apparatus. Wear protective clothing.

**SECTION 6: Accidental release measures****6.1 Personal precautions, protective equipment and emergency procedures****For non-emergency personnel**

Refer to protective measures listed in sections 7 and 8. Ensure adequate ventilation. Avoid contact with skin, eyes and clothing.

**For emergency responders**

No data available. Personal protective equipment (PPE) - see Section 8.

**6.2 Environmental precautions**

Do not discharge into the drains/surface waters/groundwater. Do not discharge into the subsoil/soil.

**6.3 Methods and material for containment and cleaning up**

Take up with absorbent material (e.g., sand, kieselguhr, universal binder). When collected, handle material as described under the section heading "Disposal considerations".

**6.4 Reference to other sections**

No data available.

**SECTION 7: Handling and storage****7.1 Precautions for safe handling****Advice on safe handling**

Risks inherent to handling the product must be minimised by applying the appropriate protective and preventive measures. Working processes should - so far as possible, according to the state of the art - be designed to rule out bodily contact or the release of hazardous substances.

**General protective and hygiene measures**

**Trade name:** CELLCLEAN**Current version :** 3.1.1, issued: 04.05.2022**Replaced version:** 3.1.0, issued: 01.06.2021**Region:** GB

Do not eat, drink or smoke during work time. Keep away from food, drink and animal feeding stuffs. Wash hands and skin before breaks and after work. Do not inhale vapours. Avoid contact with eyes and skin. Remove soiled or soaked clothing immediately. Provide eye wash fountain in work area. Have emergency shower available.

**Advice on protection against fire and explosion**

No special measures necessary.

**7.2 Conditions for safe storage, including any incompatibilities****Technical measures and storage conditions**

Keep container tightly closed in a cool, well-ventilated place. Protect from light.

**Recommended storage temperature**

Value 1 - 30 °C

**Requirements for storage rooms and vessels**

Containers which are opened must be carefully closed and kept upright to prevent leakage.

**Incompatible products**

Do not store together with: Acids; Ammonia; Reducing agents

**7.3 Specific end use(s)**

No data available.

**SECTION 8: Exposure controls/personal protection****8.1 Control parameters**

No parameters available for monitoring.

**8.2 Exposure controls****Appropriate engineering controls**

No data available.

**Personal protective equipment****Respiratory protection**

If workplace exposure limits are exceeded, a respiration protection approved for this particular job must be worn. In case of aerosol and mist formation, take appropriate measures for breathing protection in the event workplace threshold values are not specified.

**Eye / face protection**

Safety glasses with side protection shield (EN 166); Face shield

**Hand protection**

Sufficient protection is given wearing suitable protective gloves checked according to i.e. EN 374, in the event of risk of skin contact with the product. Before use, the protective gloves should be tested in any case for its specific work-station suitability (i.e. mechanical resistance, product compatibility and antistatic properties). Adhere to the manufacturer's instructions and information relating to the use, storage, care and replacement of protective gloves. Protective gloves shall be replaced immediately when physically damaged or worn. Design operations thus to avoid permanent use of protective gloves.

**Other**

Normal chemical work clothing.

**Environmental exposure controls**

No data available.

**SECTION 9: Physical and chemical properties****9.1 Information on basic physical and chemical properties**

<b>State of aggregation</b>
liquid
<b>Form/Colour</b>
liquid
chartreuse
<b>Odour</b>
slightly pungent

Trade name: CELLCLEAN

Current version : 3.1.1, issued: 04.05.2022

Replaced version: 3.1.0, issued: 01.06.2021

Region: GB

<b>pH value</b>				
Value	12.41			
<b>Boiling point / boiling range</b>				
Value	100 °C			
<b>Melting point/freezing point</b>				
No data available				
<b>Decomposition temperature</b>				
No data available				
<b>Flash point</b>				
No data available				
<b>Ignition temperature</b>				
No data available				
<b>Flammability</b>				
No data available				
<b>Lower explosion limit</b>				
No data available				
<b>Upper explosion limit</b>				
No data available				
<b>Vapour pressure</b>				
No data available				
<b>Relative vapour density</b>				
No data available				
<b>Relative density</b>				
No data available				
<b>Density</b>				
Value	1.1	-	1.2	g/ml
Reference temperature		20		°C
<b>Solubility in water</b>				
Reference temperature		20		°C
Comments	soluble			
<b>Solubility</b>				
No data available				
<b>Partition coefficient n-octanol/water (log value)</b>				
No data available				
<b>Kinematic viscosity</b>				
No data available				
<b>Particle characteristics</b>				
No data available				

**9.2 Other information**

<b>Other information</b>
No data available.

**SECTION 10: Stability and reactivity****10.1 Reactivity**

Dangerous reactions are not expected if the product is handled according to its intended use.

**10.2 Chemical stability**

Stable under recommended storage and handling conditions (See section 7).

**10.3 Possibility of hazardous reactions**

Trade name: CELLCLEAN

Current version : 3.1.1, issued: 04.05.2022

Replaced version: 3.1.0, issued: 01.06.2021

Region: GB

None, when used as directed.

**10.4 Conditions to avoid**

None, if handled according to intended use.

**10.5 Incompatible materials**

Product reacts with: Acids; Ammonia; Reducing agents

**10.6 Hazardous decomposition products**

Chlorine

**SECTION 11: Toxicological information****11.1 Information on hazard classes as defined in Regulation (EC) No 1272/2008****Acute oral toxicity**

No data available

**Acute dermal toxicity**

No data available

**Acute inhalational toxicity**

No data available

**Skin corrosion/irritation**

No data available

**Serious eye damage/irritation**

No data available

**Respiratory or skin sensitisation**

No data available

**Germ cell mutagenicity**

No data available

**Reproduction toxicity**

No data available

**Carcinogenicity**

No data available

**STOT - single exposure**

No data available

**STOT - repeated exposure**

No data available

**Aspiration hazard**

No data available

**Delayed and immediate effects as well as chronic effects from short and long-term exposure**

Inhalation of the vapours causes irritation of the respiratory tract and mucous membrane, headaches, nausea, giddiness, vomiting. After accidental ingestion, risk of oesophageal and stomach perforation! Frequent persistent contact with the skin can cause skin irritation. Eye contact with the product may lead to irritation. Influence of the product with the eyes can lead to blindness.

**11.2 Information on other hazards****Endocrine disrupting properties**

No data available.

**Other information**

No data available.

**SECTION 12: Ecological information****12.1 Toxicity****Toxicity to fish (acute)**

**Trade name:** CELLCLEAN**Current version :** 3.1.1, issued: 04.05.2022**Replaced version:** 3.1.0, issued: 01.06.2021**Region:** GB

No data available

**Toxicity to fish (chronic)**

No data available

**Toxicity to Daphnia (acute)**

No data available

**Toxicity to Daphnia (chronic)**

No data available

**Toxicity to algae (acute)**

No data available

**Toxicity to algae (chronic)**

No data available

**Bacteria toxicity**

No data available

**12.2 Persistence and degradability**

No data available.

**12.3 Bioaccumulative potential**

No data available.

**12.4 Mobility in soil**

No data available.

**12.5 Results of PBT and vPvB assessment****Results of PBT and vPvB assessment**

PBT assessment

No data available.

vPvB assessment

No data available.

**12.6 Endocrine disrupting properties**

No data available.

**12.7 Other adverse effects**

No data available.

**12.8 Other information****Other information**

Do not discharge product unmonitored into the environment.

**SECTION 13: Disposal considerations****13.1 Waste treatment methods****Product**

Allocation of a waste code number, according to the European Waste Catalogue, should be carried out in agreement with the regional waste disposal company.

**Packaging**

Residues must be removed from packaging and when emptied completely disposed of in accordance with the regulations for waste removal. Incompletely emptied packaging must be disposed of in the form of disposal specified by the regional disposer.

**SECTION 14: Transport information****14.1 Transport ADR/RID/ADN**

Class	8
Classification code	C9
Packing group	II
Hazard identification no.	80
UN number	UN1791
Proper shipping name	HYPOCHLORITE SOLUTION
Tunnel restriction code	E
Label	8

**Trade name:** CELLCLEAN**Current version :** 3.1.1, issued: 04.05.2022**Replaced version:** 3.1.0, issued: 01.06.2021**Region:** GBEnvironmentally hazardous  
substance mark

Symbol "fish and tree"

**14.2 Transport IMDG**

Class	8
Packing group	II
UN number	UN1791
Proper shipping name	HYPOCHLORITE SOLUTION
EmS	F-A, S-B
Marine Pollutant (according Index IMDG Code)	P
Label	8
Marine pollutant mark	Symbol "fish and tree"

**14.3 Transport ICAO-TI / IATA**

Class	8
Packing group	II
UN number	UN1791
Proper shipping name	Hypochlorite solution
Label	8

**14.4 Other information**

No data available.

**14.5 Environmental hazards**

Information on environmental hazards, if relevant, please see 14.1 - 14.3.

**14.6 Special precautions for user**

No data available.

**14.7 Maritime transport in bulk according to IMO instruments**

Not relevant

**SECTION 15: Regulatory information****15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture****EU regulations****Regulation (EC) No 1907/2006 (REACH) Annex XIV (List of substances subject to authorisation)**

According to the data available and/or specifications supplied by upstream suppliers, this product does not contain any substances considered as substances requiring authorisation as listed on Annex XIV of the REACH regulation (EC) 1907/2006.

**REACH candidate list of substances of very high concern (SVHC) for authorisation**

According to available data and the information provided by preliminary suppliers, the product does not contain substances that are considered substances meeting the criteria for inclusion in annex XIV (List of Substances Subject to Authorisation) as laid down in Article 57 and article 59 of REACH (EC) 1907/2006.

**Regulation (EC) No 1907/2006 (REACH) Annex XVII: RESTRICTIONS ON THE MANUFACTURE, PLACING ON THE MARKET AND USE OF CERTAIN DANGEROUS SUBSTANCES, MIXTURES AND ARTICLES**

The product is considered being subject to REACH regulation (EC) 1907/2006 annex XVII.	No 3
--	------

The product contains following substance(s) that are considered being subject to REACH regulation (EC) 1907/2006 annex XVII.

No	Substance name	CAS no.	EC no.	No
1	sodium hydroxide	1310-73-2	215-185-5	75
2	sodium hypochlorite, solution	7681-52-9	231-668-3	75

**Directive 2012/18/EU on the control of major-accident hazards involving dangerous substances**

This product is subject to Part I of Annex I, risk category:	E1
--	----

**15.2 Chemical safety assessment**

No data available.

**SECTION 16: Other information**



**Trade name:** CELLCLEAN

**Current version :** 3.1.1, issued: 04.05.2022

**Replaced version:** 3.1.0, issued: 01.06.2021

**Region:** GB

**Sources of key data used to compile the data sheet:**

Regulation (EC) No 1907/2006 (REACH), 1272/2008 (CLP) as amended in each case.

Directives 2000/39/EC, 2006/15/EC, 2009/161/EU, (EU) 2017/164.

National Threshold Limit Values of the corresponding countries as amended in each case.

Transport regulations according to ADR, RID, IMDG, IATA as amended in each case.

The data sources used to determine physical, toxic and ecotoxic data, are indicated directly in the corresponding section.

**Full text of the H- and EUH- phrases drawn up in sections 2 and 3 (provided not already drawn up in these sections)**

H318	Causes serious eye damage.
H400	Very toxic to aquatic life.
H411	Toxic to aquatic life with long lasting effects.

**Notes relating to the identification, classification and labelling of substances and mixtures ((EC) No 1272/2008, Annex VI)**

B	Some substances (acids, bases, etc.) are placed on the market in aqueous solutions at various concentrations and, therefore, these solutions require different classification and labelling since the hazards vary at different concentrations. In Part 3 entries with Note B have a general designation of the following type: 'nitric acid ... %'. In this case the supplier must state the percentage concentration of the solution on the label. Unless otherwise stated, it is assumed that the percentage concentration is calculated on a weight/weight basis.
---	---

**Creation of the safety data sheet**

UMCO GmbH

Georg-Wilhelm-Str. 187, D-21107 Hamburg

Tel.: +49 40 / 555 546 300 Fax: +49 40 / 555 546 357 e-mail: [umco@umco.de](mailto:umco@umco.de)

This information is based on our present knowledge and experience.

The safety data sheet describes products with a view to safety requirements.

It does not however, constitute a guarantee for any specific product properties and shall not establish a legally valid contractual relationship.

Document protected by copyright. Alterations or reproductions require the express written permission of UMCO GmbH.

Prod-ID 24695